

Scientific Imaging with ImageJ

Jérôme Boulanger

June 23rd, 2022

Introduction

The aim of this session is to

- improve handling of scientific images for quantification and illustration,
- understand core concepts about images,
- discover ImageJ.

Online resources are available here:

- <http://imagej.nih.gov>
- <https://imagej.net>
- <https://forum.image.sc>

Context

Quantification

- Measure objects properties (intensity, area, shape, ...)
- Count objects
- Relationship (co-occurrence, co-localisation), hierarchical.
- Speed, motion, ...

Illustration

- Use ImageJ to prepare the elements of a figure and import them in Illustrator.
- Relationship between markers
- Localisation of organelles
- Scale of objects

Image integrity

When preparing the figure:

- Apply adjustments to the entire image
- Track the sequence of manipulations (scripts/macro)

Include in the legend or methods section the following information:

- Equipment (microscopes/objective lenses, cameras, detectors, filters) and acquisition software used. microscope controlled by Zen black and equipped with a 63x/1.4 oil objective using a 488nm laser and GaAsp detectors.
- Time and space sampling (pixel size), image bit depth, temperature, imaging medium, fluorochromes
- Mention the look up table applied.
- Processing software and manipulations (deconvolution, 3D reconstructions, volume rendering, filtering, nonlinear operation, thresholding and projection).

Undocumented image manipulations can lead to accusations of research misconduct


Nature, ed. 2022. URL:

<https://www.nature.com/nature-portfolio/editorial-policies/image-integrity#microscopy>.

Office of Research Integrity, ed. 2022. URL:

<https://ori.hhs.gov/education/products/RIandImages/default.html>.

The ImageJ software



- ImageJ is a Java based image processing software.
- Java programs run on the java virtual machine (JVM)
- Can run on many operating system (Microsoft Windows, Mac OS, Linux, ...)
- Developed by Wayne Rasband since 1997 at the National Institute of Health
- ImageJ eco-system
 - ImageJ2: rewrite of ImageJ for multi dimensional data
 - Fiji: image processing package built around ImageJ2

CT Rueden et al. *BMC Bioinformatics* 18.1 (2017).





CA Schneider, WS Rasband, and KW Eliceiri. *Nature Methods* (July 2012).

Installation & Updating

First Installation

- Download Fiji from <https://imagej.net/software/fiji/downloads>.
- Unzip and save somewhere on your hard drive where default users have access.
- In Windows open the folder and double click on ImageJ-win64 and create a shortcut.

Updating

- To update ImageJ to the latest version  .
- To install/update plugins collections and manage update sites  

Fiji.app folder content

The Fiji.app folder is organized into several subfolders:

- **jars** : contains the main jar (eg ij-1.53q.jar) and extra java dependencies
- **macros** : contains macros you installed and the `StartupMacros.fiji.ijm`
- **plugins**: contains the jar (java artifacts) of the plugins
- **scripts**: contains a few matlab scripts
- **lut** : look up tables (mapping for intensity to displayed colors)

Image is data

- Sensors convert the number of detected photo-electrons into an electric voltage
- This voltage is then digitized into a number by the A/D converter.
- The image can be seen as an array of values with columns x and rows y starting at the top left corner at $(0,0)$.

Try it yourself

- Open a sample image : `File >> Open Sample >> Blobs`
- Click on the `>>` icon and tick `Pixel Inspector`




Image is data

- Sensors converts the number of detected photo-electrons into an electric voltage
- This voltage is then digitized into a number by the A/D converted.
- The image can be seen as an array of values with columns x and rows y starting at the top left corner at $(0,0)$.

Try it yourself

- Open a sample image : File Open Sample Blobs
- Click on the >> icon and tick Pixel Inspector

250	250	253	253	253	250	253	253	253	253	253	253	253	253	253	253	253	253	255	199	133
253	250	253	250	253	253	253	250	253	253	253	253	253	253	253	253	253	253	255	246	184
250	253	250	253	253	250	253	253	253	253	250	253	253	253	253	253	253	253	255	246	157
250	250	250	253	250	253	253	253	253	253	253	253	253	253	253	253	255	246	129	103	112
250	250	253	253	253	253	253	253	253	253	253	253	253	253	253	255	234	92	82	95	90
253	250	250	250	253	250	253	253	253	253	253	253	253	253	255	246	157	86	82	82	78
253	250	250	253	253	253	253	253	253	253	253	253	255	255	157	67	60	66	82	82	72
250	253	250	253	253	253	253	253	253	253	253	255	255	201	92	66	60	66	82	90	82
253	253	253	253	253	253	253	253	253	255	157	67	82	82	72	72	78	90	90	90	96
253	253	253	253	253	253	253	253	255	157	60	78	82	82	82	78	85	96	103	85	78
253	253	253	253	253	253	253	246	129	60	67	72	72	72	73	85	96	96	96	85	78
250	250	253	253	253	253	255	86	41	53	53	53	60	66	85	96	110	103	91	91	90
253	253	253	253	253	253	242	195	144	53	57	60	60	66	90	103	103	99	96	96	96
250	250	253	253	253	253	253	253	255	105	67	72	82	96	103	110	103	99	91	103	110
250	253	250	253	253	253	253	255	246	105	98	112	121	121	121	112	121	91	103	110	121
250	250	250	253	253	253	253	253	255	220	129	126	126	126	130	121	117	91	99	106	121
250	253	250	253	253	253	253	253	255	157	121	126	126	130	121	110	106	99	110	127	137
250	253	250	253	253	253	253	253	255	202	129	112	110	103	103	99	106	110	130	140	146
253	250	253	253	253	253	253	255	255	255	234	129	112	99	91	99	99	106	110	127	140
250	253	250	253	253	253	253	253	253	255	246	129	110	99	110	103	91	106	110	116	106
250	253	250	253	253	253	253	253	255	255	202	112	99	99	110	91	99	106	106	112	104
253	250	253	250	253	253	253	253	255	234	129	117	91	91	99	84	91	106	106	104	104
250	253	250	253	253	253	253	246	216	129	106	99	78	85	91	91	96	96	96	104	96

General image formats

- Digital images can be saved from the volatile memory (RAM) of a computer to its persistent memory (HDD/SSD) as files with various formats.
- Most file formats include some sort of compression:
 - Lossless compression (PNG): original values can be exactly retrieved
 - Lossy compression (JPEG): information is lost when storing the image
 - Both (TIF): some formats are containers which can include various type of compression
- ImageJ can read and write natively a few format such as TIFF, PNG, JPG
- Use TIF for saving intermediate results (with LZW compression) and PNG for figures



Microscope vendor image formats

Each microscopy company has developed a format that can store multiple series with each a multi-channel multi-plane image stack.

- Zeiss
 - LSM : TIF based file format with additional metadata and LZW compression
 - CZI : JPEG-XR and Zstd lossless compressed image
- Nikon
 - ND2 : JPEG-2000 lossless compression
- Leica
 - LIF : Customized TIF file format
- Olympus
 - OIF: multi file formats with associated images in a folder
 - OIB: store multiple OIF and dependent image in one file

Some standardized image format have also been developed to increase inter-operability:

- OME-TIFF, OME-ZARR, ICS, HDF5

Metadata

- Metadata are essential to the interpretation and processing of the image
- Acquisition software collect additional information that are stored in the files
 - Pixel size
 - Spacing between z planes
 - Detectors
 - Objective
 - Emission wavelength
 - ...
- The Open Microscopy Environment (OME) defines a specification for storing data on biological imaging.
- Findable Accessibility, Interoperability, and Reusability (FAIR) principle.

URL: <https://docs.openmicroscopy.org/ome-model/5.6.1/developers/model-overview.htm>.
Mark D. Wilkinson et al. *Scientific Data* 3 (Mar. 15, 2016). DOI: 10.1038/sdata.2016.18.
linkert`metadata`2010.

Loading an image


- Fiji will detect supported file format when selecting **File** » **Open** and use bioformat if needed to load an image from disk to memory.
- You can also use **Plugins** » **Bio-Formats** » **Bio-Formats Importer** to use Bio-formats directly
- The Bio-Formats plugin remember previous settings, so use the tool directly to define the way the drag and drop or **File** » **Open** behaves on non-native formats.

Pixel size calibration & scale bar

- Use **Image** » **Properties...** to check and set the pixel size in each dimensions
- Use **Analyze** » **Set Scale...** to use an known distance to set the scale
- Use **Analyze** » **Tools** » **Scale Bar** to display the scale on a calibrated image



Quantization

- The intensity values are quantized into grey levels
- Images are saved as 8-bit or 16-bit images by acquisition software.
 - For camera sensor (CCD, CMOS), photo-electron are accumulated in wells and converted into analogue voltages that are then digitized. Well depth is in the order of 40000 electrons.
 - For photo-multiplier tubes (PMT) photo-electrons are amplified using a chain of dynodes enabling to count individual photons.
- Use  to check & convert an image to 8-bit, using the current dynamic range.
- 32-bit mode will allow to preserve the information if the intensities get out of the initial range.

Bit depth	Dynamic range
8	0-255
12	0-4095
14	0-16383
16	0-65535

Brightness & Contrast adjustment

- **Image** » **Adjust** » **Brightness & Contrast...** or **⌘** + **↑** + **C** launches the B&C tool
- Pixels intensities are mapped linearly to the displayed intensity.
- Enable to discard unused part of the dynamic range
- Adjust the Minimum and Maximum of the dynamic range which will be displayed.
- Adjust the Brightness and Contrast to change how they are display.
- Press **Auto** to stretch the intensities between a predefined percentage of saturated pixels. Equivalent to **Process** » **Enhance Contrast...** with the normalized option left unticked.
- Press **Apply** to change the pixel values.
- When changing the image mode (16-bit to 8-bit or 32-bit to 8-bit) the B&C settings are applied to the image.

Brightness & Contrast adjustment

Ty it yourself


- Open the boat example image **File** > **Open Samples** > **Boats**
- Press the short cut **⌘** + **↑** + **C**
- Change the minimum and maximum sliders and observe the image
- Press **Set** and set the Minimum displayed value to 50 and the maximum displayed value to 200, and press the **OK** button.
- Finally press **Apply** and observe the change of the histogram in the B&C window.

Histogram

- Shows the frequency of each intensity values
- Use **Analyze** > **Histogram** or **ctrl** + **h** to display the histogram
- Helps identify clipped values (under/over exposed, saturated)
- Visualise classes of pixels

Look up table (LUT)


- Intensities are mapped to a color defined by a red, green and blue triplet value for display.
- Converting an image to RGB applies the current LUT (**Apply LUT**, apply the B&C settings)
- LUT can be inverted (black becomes white, ...)
- Use LUT to help people with colour blindness



Look up table (LUT)

Try it yourself

- Open a grey scale image `File >> Open Samples >> Boats`
- Apply a LUT to the image `Image >> Look Up Tables >> 5 Ramps`
- Convert the image to RGB `Image >> Type >> RGB Color`
- Observe the values at each pixel
- Open the blobs example `File >> Open Samples >> Blobs` and comment on the LUT
- Explore other commands `Image >> Color >> Display LUTs`, `Image >> Color >> Edit LUT`
- Open the Fluorescent Cells sample and change the red channel to magenta.



Multi-channel images

- Multi-channel images can be easily acquired by microscopes (multiple wavelength/fluorophores, brightfield/DIC, ...)
- There are two types of color images:
 - RGB images store color information into a single (8×3) 24 bit image
 - Composite images are images whose channels are stored as individual (8,16-32-bit) image planes
- On a composite image use `Image > Color > Channel Tool` or `ctrl + ↑ + z` to select the channel to display
- Create a multi channel image using `Image > Color > Merge Channels...`
- Split channels of an image using `Image > Color > Split Channels`

Try it yourself

- Open a composite image `File > Open Samples > Fluorescent Cells`
- Use the channel tool to display each colour individually

Multi-dimensional images

Stack & Hyperstacks

- Images can have many dimensions (axis) associated to them channel c, time t, depth.
- Image stacks represent up to 4 dimensions (xyzc or xytc) with a maximum of 3 channels.
- Hyperstacks can have up to 5 dimensions xyzct with no limits on the number of channels.
- Virtual stacks enable to load only the images planes which are visualised.

To quickly visualize multi-dimensional data, use :

- Orthogonal views (Image > Stack > Orthogonal Views) or **ctrl** + **↑** + **H** can help visualize 3D data
- Maximum intensity projection Image > Stack > Z Project. . .
- Reslice a stack with Image > Stack > Reslice. . . or press **/** each pixel. Activate the BIG-EPFL update site and search for “Extended Depth Of Field”

To synchronize two stacks use Analyze > Tools > Synchronize Windows

Basic manipulation

- Duplicate **ctrl** + **↑** + **D** allows to crop, select channels and slices
- Use **+** and **-** or the magnifying glass tool to zoom in and out the displayed image.
- Use the hand tool to pan within a zoomed image.
- Changing the actual pixel size:
 - **Image** » **Scale** scale the image, if “create a new image” is not ticked, the image keep the same number of pixels and is cropped or padded as necessary.
 - **Image** » **Adjust** » **Size...** scale the image in place without cropping or creating a new image.
- Use the search bar to look for tools.

Selection / Region Of Interest (ROI)

- Use the rectangle, circle, polygon, freehand, line tools to create selections.
- Use **Edit** » **Selection** » **Specify** to enter manually the coordinate of a shape
- Transfer selection across images using **Edit** » **Selection** » **Restore Selection** or **ctrl** + **↑** + **E**
- Store selection into the ROI manager **Edit** » **Selection** » **Add to manager** or **t**
- Use **Edit** » **Selection** » **Select None** **ctrl** + **↑** + **A** to be sure no selection is active
- ROI from the ROI Manager can be stored as individual proprietary .roi file or as a collection into a .zip file.
- ROIs can be combined together using AND, OR, XOR logical operation
- Use **Edit** » **Draw** and **Edit** » **Fill** to set the pixel values on the contours resp. the inside of the shape to the value define by the colour picker tool.

Threshold

- Thresholding an image gives a binary image: 8-bit image with values either 0 or 255 whether the pixel intensity lies within two values (upper & lower).
- Simplest form of image segmentation.
- Use the sort cut `ctrl` + `↑` + `T` to display the threshold tool.
 - Use the sliders to set the upper and lower values
 - Auto thresholds: Fiji includes a set of automatic threshold ¹
 - Stack histogram: the statistics will be computed on the all stack
- Use `Apply` or `Process` `»` `Binary` `»` `Convert to mask` to apply the threshold and make the image binary.
- The menu `Process` `»` `Binary` `»` `Make binary` will show a dialog if a threshold has been set.
- The LUT of binary images is inverted

¹<https://imagej.net/plugins/auto-threshold>

From thresholds to regions of interests


- Use **Edit** » **Selection** » **Create selection** to convert the threshold to a selection (can be a composite selection).
- Use the **Analyze** » **Analyze Particles...** to convert connected components into individual ROIs and add them to the ROI Manager.

Try it yourself

1. Open the blob sample image
2. Threshold the image
3. Add connected components to the ROI Manager and get the number of blobs.

Measurements


- Use **Analyze >> Set Measurements** to define the quantity to be measured
- Select an individual ROI and press **ctrl + M** to add measurements to the Result table.
- To measure the properties of all ROI stored in the ROI Manager, select them all with **ctrl + A** and press **Measure**



Binary image processing

Binary image processing relies mostly on “morphological image processing” which process image by shapes. It can help refine the mask generated by thresholding:

- **Process > Binary > Dilate** grow the masks, similar to **Process > Filter > Maximum...**
- **Process > Binary > Erode** shrink the masks, similar to **Process > Filter > Minimum...**
- **Process > Binary > Open** is equivalent to erode followed by dilate, removes small objects and lines
- **Process > Binary > Close-** is equivalent to dilate followed by erode, bridge non touching objects etc
- **Process > Binary > Watershed** separate touching objects
- **Process > Binary > Fill Holes** will fill holes in the masks



Filtering

Image filtering can help reduce noise


- Box filtering `ctrl` + `↑` + `S`
- Gaussian filtering `Process` `Filters` `Gaussian Blur...`
- Median filtering `Process` `Filters` `Gaussian Blur...`

Try it yourself

- Create an image with a light grey square on a darker grey square.
- Add noise to the image
- Launch the histogram and press `Live`
- Apply some filtering to the image
- Apply a threshold

Background correction

- Background level can vary for several reason
 - Uneven illumination
 - Scattering
 - Auto fluorescence
- Simple generic correction
 - Rolling ball from Process Subtract Background...
 - Top Hat (image - grayscale opening) from Process Filter Top Hat...
- Illumination correction
 - BASIC plugin




Drift correction

- Install Stackreg from the BIG-EPFL update site
- Plugin » Registration » Stackreg
- The plugin can register drift using several motion models such as
 - Translation
 - Rigid body
 - Rotation
 - Scaled rotation
 - Affine

Line profile


- Line profile enable to measure distances visually.
- Draw a line and select **Analyze** > **Plot Profile** or press **ctrl** + **K**
- Use **Data** > **Add fit** to approximate the data with a model



Kymograph


1. Create a maximum intensity projection
2. Select lines using a freehand tool and add them to the ROI Manager
3. Double click on the line tool to access the line width dialog
4. Reslice the original stack

Velocities can be then extracted from the kymographs using a macro for example.



Pixel classifier

- Weka is a machine learning library in Java
- The weka plugin allows to annotate, train, export and apply the classifier




Nuclei segmentation with StarDist

- StartDist is deep learning based approach based on star convex polygons
- It has been trained to identify cell nuclei.
- Install with **Help** > **Update** and locate the StarDist and CSBDeep entries
- Launch trackmate from **Plugins** > **StarDist**



Tracking with TrackMate

- Detect spots and link them overtime
- Install trackmate using **Help** > **Update** and locate the TrackMate entry
- Launch trackmate from **Plugins** > **Tracking** > **Trackmate**
- Trackmate can also use stardist or weka for detecting objects



Basic quantification workflow

1. Open the image
2. Filtering
3. Threshold
4. Refine binary image
5. Extract ROIs
6. Measure

Basic figure preparation workflow

1. Open the image
2. Maximum intensity projection
3. Adjust contrasts
4. Convert to RGB
5. Add a scale bar
6. Export as PNG

Figure preparation tools

Manually

- Convert to RGB, add a scale bar without text, import in Illustrator

figureJ

- Update sites IBMP-CNRS, ImageScience

EZFig

- Update sites: EZF

QuickFigure

- Update sites: QuickFigure

Extension

ImageJ functionality can be extended and customised using:

- Macros (.txt, .ijm files)
- Scripts (Javascript, Python, Groovy . . .)
- Plugins (in Java)

ImageJ Cheatsheet

Image properties

ctrl + ↑ + P

Duplicate/crop

ctrl + ↑ + D

BrightnessContrast

⌘ + ↑ + C

Channel Tools

ctrl + ↑ + Z

Add to ROI Manager

t

Restore Selection

ctrl + ↑ + E

Select All

ctrl + A

Select None

ctrl + ↑ + A

Histogram

ctrl + H

Orthoslices

ctrl + ↑ + H

Reslice

/

Threshold

ctrl + ↑ + T

Measure

ctrl + M

Smooth

ctrl + ↑ + S

Plot profile

ctrl + K